



The next generation in Point Of Care monitoring...



FreeStyle Optium Neo H



The dedicated hospital blood glucose meter

The FreeStyle NeoH can perform both blood glucose and blood Ketone tests. Ensure to familiarise yourself with both modules.

Click on the module below.



FreeStyle

Optium Neo H

Blood Glucose and Ketone Monitoring System

An Introduction to FreeStyle Optium NeoH




FreeStyle

Optium Neo **H**

Blood Glucose and Ketone Monitoring System

The FreeStyle Optium Neo H System at a Glance



FreeStyle

Optium Neo H

Blood Glucose and Ketone Monitoring System

The FreeStyle Optium Neo H Blood Glucose Test Strip

The FreeStyle Optium Neo H meter is designed to be used only with FreeStyle Optium Neo H Blood Glucose Test Strips.

- 0.6µL sample size
- 5-second result countdown
- Assay range 1.1-27.8 mmol/L
- Haematocrit range 15-65%
- Top or tip sample application
- Visual confirmation of sufficient sample
- Fill trigger ensures sufficient sample
- Sample Types – Capillary, venous, arterial & neonatal
- Each pack contains 100 individually foil wrapped test strips, calibration bar and test strip insert



FreeStyle

Optium Neo H

Blood Glucose and Ketone Monitoring System

Meter calibration

- The calibration procedure programmes the meter with the lot number, expiry date and test strip technology
- Always calibrate the meter with each new pack of test strips – failure to calibrate properly will cause incorrect results
- With the lot number facing toward you, insert the contact bars of the calibrator into the meter. The meter turns on automatically
- LOT and the lot number appear in the display window
- Check that the lot number on the meter display window matches with the number on the test strip calibrator and the last five digits on both the test strip foil packet and test strip insert
- Retain the calibration bar for the duration of that pack
- Calibration is complete
- Only use the calibrator that is packaged in the pack of test strips you are using



FreeStyle

Optium Neo H

Blood Glucose and Ketone Monitoring System



FreeStyle

Optium Neo H

Blood Glucose and Ketone Monitoring System

Performing a Patient Test

The FreeStyle Optium Neo H system enables easy and accurate bedside testing of blood glucose levels

Blood glucose testing

- ▶ 5 second test time
- ▶ Small 0.6 μ L sample size
- ▶ Calibrator provides expiry check and calibrates the meter
- ▶ Control solution mark

Patient Testing

1. Prepare the sample site:
 - i. Before you obtain a blood drop, make sure the sample site is clean, dry and warm
 - ii. Hang the arm down before lancing the sample site, to help aid blood flow
2. Open the individually foil wrapped test strip at the notch and tear down
3. Insert the test strip into the meter until it stops. This will turn the meter on.
4. Check that the LOT number of test strip you are using matches the LOT number shown on the meter display
5. The  blinks, indicating the meter is ready for you to apply a sample to the test strip
 - i. Ensure you follow organisational policy for sample collection



FreeStyle

Optium Neo 

Blood Glucose and Ketone Monitoring System

Patient testing

6. If testing glucose on venous blood press the **V** to mark the test until  appears indicating that the meter is now ready for you to apply venous blood to the test strip
7. Apply a sample with the meter in a horizontal (flat) position to avoid it entering the test strip port
8. If you are applying the sample from a syringe, apply away from the direction of the test strip port



FreeStyle

Optium Neo H

Blood Glucose and Ketone Monitoring System

Patient testing

9. Hold sample to the test strip until 3 short lines appear on the screen



10. The results will be displayed after a 5-second countdown



FreeStyle

Optium Neo H

Blood Glucose and Ketone Monitoring System

Patient testing

When the result appears on the meter screen the test is complete.

The results will be stored in the meter's memory.

Press and hold the  button to turn off the meter.



FreeStyle

Optium Neo H

Blood Glucose and Ketone Monitoring System

Summary for blood glucose testing



1. Inserts test strip,  blinks
2. Apply blood
3. Sample detected, 5 seconds countdown begins
4. Result is displayed



FreeStyle

Optium Neo H

Blood Glucose and Ketone Monitoring System

Out of range indicators

High blood glucose results

- A result falling outside the range set by your organisation may require some intervention from the operator, such as repeating the test
 - The meter will display a solid yellow arrow
- For guidance on out of range results, please follow your organisational policy



Low blood glucose results

- A result falling outside the range set by your organisation will require some intervention from the operator, such as repeating the test
 - The meter will display a solid red arrow
- For guidance on out of range results, please follow your organisational policy



FreeStyle

Optium Neo H


Blood Glucose and Ketone Monitoring System

Display Messages



LO

Result is less than
1.1 mmol/L



HI

Glucose result is
above 27.8 mmol/L



KET

Result is higher than or
equal to 13.3 mmol/L



FreeStyle

Optium Neo **H**

Blood Glucose and Ketone Monitoring System

Results & Memory

The test is complete when the result appears on the meter screen.

- The result is stored in the meter's logbook
- The meter logbook can store up to 1,000 events – including blood glucose, blood ketone, and control solution results, and other meter information
- The logbook can be accessed by going to the Home Screen and selecting 
- To recall patient results:

With the FreeStyle Optium Neo H turned on:

1. Press the  to view logbook events
2. Press  or  to scroll through log book events



While on the Home screen, press  to view logbook events




FreeStyle

Optium Neo H


Blood Glucose and Ketone Monitoring System

Logbook: View Logbook Events


Press  to scroll through to view logbook events



Blood Glucose Result



Venous Blood Result



Blood Ketone Result


Control solution test results will alternate between 'CtL' and result



Glucose Control
Solution Result



Ketone Control
Solution Result



FreeStyle

Optium Neo H

Blood Glucose and Ketone Monitoring System

Logbook: View Blood Glucose Averages

1. While in the logbook, press  any time to view the 7-day average

2. Press  or  to scroll through 7-, 14- and 30-day averages. Press  to return to logbook events.



FreeStyle

Optium Neo H

Blood Glucose and Ketone Monitoring System

Error Messages

Message	What it means
E1	Temperature is too hot or too cold for the meter to work correctly
E2/E9	Meter error
E3	<ul style="list-style-type: none">i. Blood drop is too smallii. Incorrect test procedureiii. There may be a problem with the test strip
E4	<ul style="list-style-type: none">i. The blood glucose level may be too high to be read by the systemii. There may be a problem with the test strip
E5	<ul style="list-style-type: none">i. Blood was applied to the test strip too soon
E6	<ul style="list-style-type: none">i. Calibration errorii. Test strip error
E7	<ul style="list-style-type: none">i. Test strip may be damaged, used or the meter does not recognise it



FreeStyle

Optium Neo **H**

Blood Glucose and Ketone Monitoring System

Performing a Quality Control Test



FreeStyle

Optium Neo H

Blood Glucose and Ketone Monitoring System

Performing a Quality Control Test

- If QC testing is NOT performed, misleading results may lead to inappropriate treatment
- Always follow your organisational policy to avoid putting patients at risk
- The FreeStyle Optium Neo H meter may prompt when this is due
- Only use MediSense Glucose and Ketone Control Solutions and keep the meter flat to prevent excess solution running into the meter
- QC Solution is stable for 90 days following opening or until the expiry date printed on the solution bottles, whichever is sooner. Discard and do NOT use control solution past its expiry date
- Check results against the test strip insert inside the pack of test strips being used as QC ranges will vary between test strip LOTs
- It is the responsibility of all users to ensure the meter has been tested before patient use



FreeStyle

Optium Neo H

Blood Glucose and Ketone Monitoring System



Performing Quality Control (QC) Tests

1. Check that the LOT number printed on the test strip and instructions for use match
2. Open the individual foil test strip at the notch and tear down
3. Insert the test strip into the meter which will turn the meter on
4. When the  blinks press the  button until this icon appears 
5. The meter is now ready for you to apply the low level control solution to the test strip
 - i. Invert the bottle to mix the solution
 - ii. Wipe the QC nozzle with clean gauze before and after each test
 - iii. Replace the cap securely on the bottle immediately after use
7. Apply a drop of the control solution with the meter in a horizontal position to avoid it entering the test strip port



Performing Quality Control (QC) Tests

8. Apply a drop of the control solution to the test strip until 3 short lines appear on the screen
9. The results will be displayed after a 5-second countdown
 - i. Compare the control solution results to the range printed on the blood glucose test strip insert; the results should fall within this range
 - i. The results will be stored in memory as a control solution result
 - ii. Record the result in the QC record book
10. If results are out-of-range
 - i. Repeat the test
 - ii. If the repeated result is out-of-range follow organisational policy for out-of-range results
11. Remove the test strip, repeat the procedure for high level control solution
12. Record HI and LO QC test result in the QC record book



FreeStyle

Optium Neo H

Blood Glucose and Ketone Monitoring System

Quality Control Testing Summary



1. Insert test strip,  blinks.
2. Press and hold the  until  appears.
3. The meter is now ready for the operator to apply control solution to the test strip.
4. Apply control solution.
5. Sample detected, countdown begins.
5. Result is displayed.

FreeStyle

Optium Neo H

Blood Glucose and Ketone Monitoring System



Maintenance



FreeStyle

Optium Neo **H**

Blood Glucose and Ketone Monitoring System

Maintenance of the meter

1. Gently clean the exterior of the FreeStyle Optium Neo H meter with a damp cloth
 - i. Acceptable cleaning solutions include mild detergent and water, 70% isopropyl alcohol or a mixture of 1 part household bleach and 9 parts water
 - ii. Do not immerse the meter in water or other liquids. Avoid getting dust, blood, control solution and any other substances in the test strip port, USB port and battery compartment
 - iii. Allow the meter to dry
2. Do not try to clean the test strip port.
3. Do not pour liquid into the test strip port or onto the buttons.
4. Do not immerse the meter in water or any other liquid.



FreeStyle

Optium Neo H

Blood Glucose and Ketone Monitoring System

Maintenance of the meter

1. Replace coin cell batteries when 'Low Battery Power' message is displayed
 1. The meter requires 2 replaceable CR2032 coin cell batteries.
 2. To change the batteries turn the meter over and slide the battery door (located on the side of the device) open, as shown.
2. Perform QC test as per organisational policy
3. Refer to troubleshooting guide in FreeStyle Optium Neo H operator's manual for error messages
4. When not in use, always store the FreeStyle Optium Neo H meter as per organisational policy



FreeStyle

Optium Neo H

Blood Glucose and Ketone Monitoring System



FreeStyle

Optium Neo H

Blood Glucose and Ketone Monitoring System

Blood β -Ketone Testing with FreeStyle Optium Neo H

The FreeStyle Optium Neo H system enables easy and accurate bedside testing of blood ketone levels

Blood ketone testing

- ▶ 10 second test time
- ▶ 1.5 μ L sample size
- ▶ Purple calibrator calibrates the meter
- ▶ Control solution mark

An Introduction to FreeStyle Optium NeoH



FreeStyle

Optium Neo **H**

Blood Glucose and Ketone Monitoring System

The FreeStyle Optium Neo H System at a Glance



FreeStyle

Optium Neo H

Blood Glucose and Ketone Monitoring System

FreeStyle Optium H Blood β -Ketone Test Strips

The FreeStyle Optium H Blood β -Ketone Test Strip:

- 1.5 μ L sample size
- 10-second result countdown
- Assay Range 0.0-8.0 mmol/L
- Haematocrit Range 30-60%
- Top or tip sample application
- Visual confirmation of sufficient sample
- Fill trigger ensures sufficient sample
- Sample Type – capillary (finger) and venous whole blood samples
- Each pack contains 10 individually foil wrapped test strips, calibration bar and test strip insert



FreeStyle

Optium Neo H

Blood Glucose and Ketone Monitoring System

β -Ketone Benefits

Clinical accuracy

- Measures β -hydroxybutyrate (the most significant ketone body)
- A quantitative vs. qualitative measure
- Results unaffected by many common medications

Patient convenience/compliance

- Blood vs. urine sampling
- Foil wrapped test strips
- Test integrated within blood glucose meters



FreeStyle

Optium Neo **H**

Blood Glucose and Ketone Monitoring System

Meter calibration

- The calibration procedure programmes the meter with the lot number, expiry date and test strip technology
- Always calibrate the meter with each new pack of test strips – failure to calibrate properly will cause incorrect results
- With the lot number facing toward you, insert the contact bars of the calibrator into the meter. The meter turns on automatically
- LOT and the lot number appear in the display window
- Check that the lot number on the meter display window matches with the number on the test strip calibrator and the last five digits on both the test strip foil packet and test strip insert
- Retain the calibration bar for the duration of that pack
- Calibration is complete
- Only use the calibrator that is packaged in the pack of test strips you are using



FreeStyle

Optium Neo H

Blood Glucose and Ketone Monitoring System

Performing a Patient Test



FreeStyle

Optium Neo H

Blood Glucose and Ketone Monitoring System

Patient Testing

1. Prepare the sample site:
 - i. Before you obtain a blood drop, make sure the sample site is clean, dry and warm
 - ii. Hang the arm down before lancing the sample site, to help aid blood flow
2. Open the individually foil wrapped test strip at the notch and tear down
3. Insert the test strip into the meter until it stops. This will turn the meter on.
4. Check that the LOT number of test strip you are using matches the LOT number shown on the meter display
5. The  blinks, indicating the meter is ready for you to apply a sample to the test strip
 - i. Ensure you follow organisational policy for sample collection



FreeStyle

Optium Neo 

Blood Glucose and Ketone Monitoring System

Patient testing

6. If testing ketones on venous blood press the **V** to mark the test until  appears indicating that the meter is now ready for you to apply venous blood to the test strip
7. Apply a sample with the meter in a horizontal (flat) position to avoid it entering the test strip port
8. If you are applying the sample from a syringe, apply away from the direction of the test strip port



FreeStyle

Optium Neo **H**

Blood Glucose and Ketone Monitoring System

Patient testing

9. Hold sample to the test strip until 3 short lines appear on the screen



10. The results will be displayed after a 10-second countdown



FreeStyle

Optium Neo H

Blood Glucose and Ketone Monitoring System

Patient testing

When the result appears on the meter screen the test is complete.

The results will be stored in the meter's memory.

Press and hold the  button to turn off the meter.



FreeStyle

Optium Neo H

Blood Glucose and Ketone Monitoring System

Display Messages

Message	Explanation
HI	Result above 8.0 mmol/L

This appears when a result is higher than 8.0 mmol/L
This suggests very high blood ketone levels or that there may be an error with the test strip
Repeat the test with a new test strip
This requires immediate action as defined by your facility's policies and procedures



FreeStyle

Optium Neo **H**

Blood Glucose and Ketone Monitoring System

Blood Ketone Testing summary



1. Insert ketone test strip  blinks and **KET** appears

2. Apply blood

3. Sample detected, 10 seconds countdown begins

4. Result is displayed



FreeStyle

Optium Neo H

Blood Glucose and Ketone Monitoring System

Error Messages

Message	What it means
E1	Temperature is too hot or too cold for the meter to work correctly
E2/E9	Meter error
E3	<ul style="list-style-type: none">i. Blood drop is too smallii. Incorrect test procedureiii. There may be a problem with the test strip
E4	<ul style="list-style-type: none">i. The blood glucose level may be too high to be read by the systemii. There may be a problem with the test strip
E5	<ul style="list-style-type: none">i. Blood was applied to the test strip too soon
E6	<ul style="list-style-type: none">i. Calibration errorii. Test strip error
E7	<ul style="list-style-type: none">i. Test strip may be damaged, used or the meter does not recognise it



FreeStyle

Optium Neo **H**

Blood Glucose and Ketone Monitoring System

Performing a Quality Control Test



FreeStyle

Optium Neo **H**

Blood Glucose and Ketone Monitoring System

Performing a Quality Control Test

- If QC testing is NOT performed, misleading results may lead to inappropriate treatment
- Always follow your organisational policy to avoid putting patients at risk
- The FreeStyle Optium Neo H meter may prompt when this is due
- Only use MediSense Glucose and Ketone Control Solutions and keep the meter flat to prevent excess solution running into the meter
- QC Solution is stable for 90 days following opening or until the expiry date printed on the solution bottles, whichever is sooner. Discard and do NOT use control solution past its expiry date
- Check results against the test strip insert inside the pack of test strips being used as QC ranges will vary between test strip LOTs
- It is the responsibility of all users to ensure the meter has been tested before patient use



FreeStyle

Optium Neo **H**

Blood Glucose and Ketone Monitoring System



Performing Quality Control (QC) Tests

1. Check that the LOT number printed on the test strip and instructions for use match
2. Open the individual foil test strip at the notch and tear down
3. Insert the test strip into the meter which will turn the meter on
4. When the  blinks press the  button until this icon appears 
5. The meter is now ready for you to apply the low level control solution to the test strip
 - i. Invert the bottle to mix the solution
 - ii. Wipe the QC nozzle with clean gauze before and after each test
 - iii. Replace the cap securely on the bottle immediately after use
7. Apply a drop of the control solution with the meter in a horizontal position to avoid it entering the test strip port



*KET will appear on the screen as a purple ketone test strip has been inserted



FreeStyle

Optium Neo H

Blood Glucose and Ketone Monitoring System



Performing Quality Control (QC) Tests

8. Apply a drop of the control solution to the test strip until 3 short lines appear on the screen
9. The results will be displayed after a 10-second countdown
 - i. Compare the control solution results to the range printed on the blood glucose test strip insert; the results should fall within this range
 - i. The results will be stored in memory as a control solution result
 - ii. Record the result in the QC record book
10. If results are out-of-range
 - i. Repeat the test
 - ii. If the repeated result is out-of-range follow organisational policy for out-of-range results
11. Remove the test strip, repeat the procedure for high level control solution
12. Record HI and LO QC test result in the QC record book



FreeStyle

Optium Neo H

Blood Glucose and Ketone Monitoring System

Quality Control Testing Summary



1. Insert test strip,  blinks.
2. Press and hold the  until  appears.
3. The meter is now ready for the operator to apply control solution to the test strip.
4. Apply control solution.
5. Sample detected, countdown begins.
5. Result is displayed.

FreeStyle

Optium Neo H

Blood Glucose and Ketone Monitoring System

Maintenance



FreeStyle

Optium Neo H

Blood Glucose and Ketone Monitoring System

Maintenance of the meter

1. Gently clean the exterior of the FreeStyle Optium Neo H meter with a damp cloth
 - i. Acceptable cleaning solutions include mild detergent and water, 70% isopropyl alcohol or a mixture of 1 part household bleach and 9 parts water
 - ii. Do not immerse the meter in water or other liquids. Avoid getting dust, blood, control solution and any other substances in the test strip port, USB port and battery compartment
 - iii. Allow the meter to dry
2. Do not try to clean the test strip port.
3. Do not pour liquid into the test strip port or onto the buttons.
4. Do not immerse the meter in water or any other liquid.



FreeStyle

Optium Neo H

Blood Glucose and Ketone Monitoring System

Maintenance of the meter

1. Replace coin cell batteries when 'Low Battery Power' message is displayed
 1. The meter requires 2 replaceable CR2032 coin cell batteries.
 2. To change the batteries turn the meter over and slide the battery door (located on the side of the device) open, as shown.
2. Perform QC test as per organisational policy
3. Refer to troubleshooting guide in FreeStyle Optium Neo H operator's manual for error messages
4. When not in use, always store the FreeStyle Optium Neo H meter as per organisational policy



FreeStyle

Optium Neo H

Blood Glucose and Ketone Monitoring System

Congratulations on completing this module



FreeStyle

Optium Neo **H**

Blood Glucose and Ketone Monitoring System

FreeStyle blood glucose meters assist in the monitoring of blood glucose (and ketone) levels for people diagnosed with diabetes mellitus. People diagnosed with diabetes with elevated blood ketone levels should seek medical advice. Information contained herein for distribution outside of the U.S. Only. FreeStyle and related brand marks are trademarks of Abbott Diabetes Care Inc. in various jurisdictions. All other trademarks are the property of their respective owners. Abbott Diabetes Care Australia, 666 Doncaster Road, Doncaster VIC 3180. For further information call Customer Service on 1800 801 478 or visit www.myfreestyle.com.au MSE160420044321