

EB-Pro

Do More | Faster | Focused

Wireless Flagship(192 crystal array/32 channels)/ IF winner 2021



Greater accuracy & earlier detection compared. Suitable for Bovine, Equine, Ovine, Swine , widely used in uterus, ovaries, fetal structure, fetal sex, reproductive tract disease, horse's far limb, tendon, bone, eyes, abdomen and chest organs.

EB-Pro

Do More | Faster | Focused

EB-Pro specifically engineered for all aspects of bovine and equine reproductive imaging, as well as for basic equine tendon imaging.

- ◆ IP54 performance
- ◆ Intelligent HD imaging : 32 channels, 128 elements
- ◆ 12 hours continuous working time, video glasses up to 4 hours (option)
- ◆ Special patented BMV's veterinary software
- ◆ Fast connection and conversion:
- ◆ Supports (Android, Windows) tablet, mobile phone, laptop
- ◆ Tablets and TVs can be connected at the same time, connecting faster and more suitable for conversion
- ◆ Plenty of accessories options: introducer stick, Binocular lenses, smart devices, backpack, belt, wrist straps, portable suitcase, mobile charger, wireless remote control, and air print etc



Technical Specifications

	EB-Pro
Size, weight	185 x 80 x 35 mm, 650g
Dust and splash proof	Yes
	Broadband linear rectal (Acoustic array length 64mm) Broadband rectal curved array transducer (for deeper imaging)
Elements	192 crystal array
Channels	32 channels
Frequency Max	2-11 MHz
Depth Imaging	31cm
Modes	B, M
Viewing device	Android Smartphone/ Tablet/I-Scan® Goggles.
Battery life (hrs)	10~12 hours
Battery Charging	Battery can be charged from DC power adapter or car 12v accessory socket