

## Oscillating Prophy Angle - Disposable



R&R Oscillating Disposable Prophy Angle, are supplied with soft rubber prophy cups in bags of 100/pk. The R&R's are pre-lubricated and iM3's unique blue colour is your guarantee of the prophy angle lasting a complete procedure! The 60° reciprocating 'back and forth' action of the prophy cup prevents muzzle hair from being tangled and allows the cup to have continuous contact on the tooth surface. The reciprocating action also reduces heat and friction on the tooth and doesn't splatter prophy paste. Guaranteed not to catch hair or whiskers!

The new L7466 Blue prophy angle supersedes the L7465 previously sold by iM3.

