

AET-V1

Outer packing box I
Size: (L)1495*(W)95*(H)155mm
Standard weight: 6.5kg

Outer packing box II
Size: (L)665*(W)585*(H)470mm
Standard weight: 27kg

Endoscope pole assembly holds up to 2 flexible endoscopes.

Monitor mount supports OR-160 and OR-100 processor or monitors, swivels 180° left and right, and offers adjustable pitch tilt.

All 4 wheels are lockable.

Endoscope pole
with a 420 mm adjustable
height range
Max Weight = 6kg

Monitor mount
Max Weight = 10kg

Top tray
455(W)*420(D)
Max Weight = 13kg

Second shelf
450(W)*400(D)
Max Weight = 10kg

Third shelf
450(W)*400(D)
Max Weight = 10kg

