

# BUSTER Easy-ID Dog Nylon Muzzle

Line of high-quality colour coded nylon muzzles that have a side easy-release snap buckle


- High quality water-resistant
- Colour coded for easy recognition
- Extra padding around the nose
- Muzzles fit securely with adjusting straps for a snug fit
- Lightweight nylon
- Reinforced stitching to ensure long-lasting use
- Available in 10 sizes from 2XS to 5XL



**COLOUR CODED  
FOR EASY SIZE IDENTIFICATION**

### Selecting the right muzzle for your dog

- Use a string to measure around your dog's nose at the position shown
- Place your index finger under the string to ensure enough space for the dog to be comfortable. Do not tighten the string; it should be free to move around the dog's nose
- Lay the string along a ruler to get appropriate length, then select the proper muzzle  
The muzzle must not be tight. It is better a little loose than too tight



### Weight/Breed Guideline

Cat. No	Size	Patient Weight Range (kg)	Patient Weight Range (lbs)	Snout Circumference (cm)	Snout Circumference (")	Colour	Breeds
279465	2XS	< 2 kg	< 5 lbs	9 cm	3½"	●	Chihuahua
279466	XS	2-4 kg	5-10 lbs	11.5 cm	4½"	●	Toy Poodle
279467	S	4-11 kg	10-25 lbs	13 cm	5¼"	●	Dachshound
279468	M	11-18 kg	25-40 lbs	15 cm	6"	●	Cocker Spaniel
279469	L	18-31 kg	40-70 lbs	16.5 cm	6½"	●	Bull Terrier
279470	XL	31-45 kg	70-100 lbs	20 cm	8"	●	Doberman
279471	2XL	45-68 kg	100-150 lbs	25 cm	10"	●	German Shepherd
279472	3XL	45-54 kg	100-120 lbs	28 cm	11"	●	Labrador Retriever
279473	4XL	54-68 kg	120-150 lbs	30 cm	12"	●	Great Dane
279475	5XL	> 68 kg	> 150 lbs	35 cm	14"	●	Giant Breeds
279477	XXS-5XL	Set of 10 BUSTER Easy-ID Dog Nylon Muzzles				Ass.	All
279487	S-2XL	Set of 5 BUSTER Easy-ID Dog Nylon Muzzles				Ass.	All
279489	XS – 2XL	Set of 6 BUSTER Easy-ID Dog Nylon Muzzles + 1 BUSTER muzzle for cat (size M)				Ass.	All

