




NOVATEN[®]

Synthetic Tendon

Veterinary tendon repair.

Characteristics

- Medical Grade UHMWPE
- Strong abrasion resistance
- Porosity enhancing tissue-in-growth
- Enabling mini-invasive surgery

Applications

For all replacement or reinforcement of small animals tendons.

Among others an operative technique has been developed for reconstruction of the Gastrocnemius Tendon.

“NOVATEN enables the reinforcement or the replacement of small animals tendons.

The main application developed is regarding the Gastrocnemius Tendon injuries in dogs but it can be used for the reconstruction of approximately any torn tendons.

”

NOVATEN[®]

Synthetic Tendon

Our implants are available in different sizes for being perfectly adapted to the bones of the animals. This enables to keep the surgery as minimally invasive and conservative as possible to preserve the proprioception.

Made of medical grade fibre (UHMWPE) which already proved their perfect biocompatibility, the implant have a stronger resistance than the original tendon.

Product Data

Designation	References	Characteristics
NOVATEN 4000	NO-4000T	4000 Newton Resistance. 4mm wide Diameter of Bone Tunnels : 3,0mm
NOVATEN 8000	NO-8000T	8000 Newton Resistance. 8mm wide Diameter of Bone Tunnels : 3,6mm

Related Products



Instrumentation

For ligament/tendon repair

A dedicated set of instrumentation has been developed for the implants NOVALIG & NOVATEN.



Interference Screw

Soft Tissue

Interference Screws have been designed for a strong mechanical fixation in ligament/tendon surgery.

CONTACT US

C/o MonacoTech
Zone F, Entrée B
4-6 avenue Albert II
98000 Monaco

Phone: +33 (0)6 07 93 77 11
Mail: info@novetech-surgery.com
Web: www.novetech-surgery.com