



IMPLANTS

End-Threaded Positive Profile ESF Pins

Micro ESF Pins - Stainless Steel



005089



End - Threads



005089



001396



001801



001398



001803

Micro ESF Pins - Stainless Steel

MWI#	SECUROS#	SHANK DIAMETER (")	THREAD DIAMETER (")	THREAD LENGTH (")	ACRYLIC THREAD LENGTH (")	OVERALL LENGTH (")
049503	005089	.045	.055	.6	1.0	3
049504	005090	.062	.071	.6	1.0	3
049505	005091	5/64	.091	.6	1.0	3
049506	005092	3/32	.114	.8	1.0	3

Mini ESF Pins - Stainless Steel

MWI#	SECUROS#	SHANK DIAMETER (")	THREAD DIAMETER (")	THREAD LENGTH (")	OVERALL LENGTH (")	RECOMMENDED DRILL BIT (mm)
028970	001396	1/16	5/64	.5	3.5	1.5

Small ESF Pins - Stainless Steel

MWI#	SECUROS#	SHANK DIAMETER (")	THREAD DIAMETER (")	THREAD LENGTH (")	OVERALL LENGTH (")	RECOMMENDED DRILL BIT (mm)
028974	001801	3/32	.132	1.0	3.5	2.4

Medium ESF Pins - Stainless Steel

MWI#	SECUROS#	SHANK DIAMETER (")	THREAD DIAMETER (")	THREAD LENGTH (")	OVERALL LENGTH (")	RECOMMENDED DRILL BIT (mm)
028972	001398	1/8	.157	1.5	5.0	3.2

Large ESF Pins - Stainless Steel


MWI#	SECUROS#	SHANK DIAMETER (")	THREAD DIAMETER (")	THREAD LENGTH (")	OVERALL LENGTH (")	RECOMMENDED DRILL BIT (mm)
028976	001803	5/32	3/16	2.3	7.0	4.0

Center-Threaded Positive Profile ESF Pins


001397

001802

001399

001804

Center - Threads

Mini ESF Pins - Stainless Steel

MWI#	SECUROS#	SHANK DIAMETER (")	THREAD DIAMETER (")	THREAD LENGTH (")	OVERALL LENGTH (")	RECOMMENDED DRILL BIT (mm)
028971	001397	1/16	5/64	1.0	5.0	1.5

Small ESF Pins - Stainless Steel

MWI#	SECUROS#	SHANK DIAMETER (")	THREAD DIAMETER (")	THREAD LENGTH (")	OVERALL LENGTH (")	RECOMMENDED DRILL BIT (mm)
028975	001802	3/32	.132	1.0	5.0	2.4

Medium ESF Pins - Stainless Steel

MWI#	SECUROS#	SHANK DIAMETER (")	THREAD DIAMETER (")	THREAD LENGTH (")	OVERALL LENGTH (")	RECOMMENDED DRILL BIT (mm)
028973	001399	1/8	.157	1.5	7.0	3.2

Large ESF Pins - Stainless Steel

MWI#	SECUROS#	SHANK DIAMETER (")	THREAD DIAMETER (")	THREAD LENGTH (")	OVERALL LENGTH (")	RECOMMENDED DRILL BIT (mm)
028977	001804	5/32	3/16	2.3	8.0	4.0