



# Aesculap Vega II

Cordless Trimmer for Vets

# Aesculap Cordless Clipper Vega II

The completely redesigned Vega II Veterinary Animal Shearer with 3 speed settings impresses with a new, powerful rechargeable battery for up to 3 hours of running time, as well as with an electronic „Speed Control“ power control that reliably shears off any type of coat, no matter how difficult.

- very good cutting results thanks to blade set made from a special steel alloy (stainless steel)
- fits the hand perfectly thanks to slim shape
- works very quietly and without vibration
- can also be operated with charging cable connected
- blade set can be removed without tools for easy cleaning
- fine-toothed precision blade set with integrated 3-stage cutting length adjustment (0.5/1.5/2.5 mm)


Speed Control: electronic speed control depending on the coat conditions for uniform shearing speed

Included in delivery: 1 x cordless clippers, 1 x clipper head (GT606), 1 x brush, 1 x charging station, 1 x Aesculap oil (GT604), 1 x instruction manual

| Technical data                | GT412   |
|-------------------------------|---|
| Blade set system              | Aesculap Special  |
| Drive                         | DC motor  |
| Stroke rate                   | approx. 5.500 strokes/min / approx. 6.000 strokes/min / approx. 6.500 strokes/min |
| Battery technology            | Li-Ion  |
| Battery capacity              | 3.200 mAh   |
| Battery run time              | 180 min   |
| Battery charge time           | 130 min   |
| Cutting width                 | 40 mm   |
| Cutting length                | 0,5 mm / 1,5 mm / 2,5 mm  |
| Speed levels                  | 3 speed levels  |
| Sound pressure level          | 60 dB   |
| Battery charge status display | 5-tier  |
| Dimensions                    | 17 cm x 4,6 cm x 4,6 cm   |
| Weight                        | 250 g   |

| Ref. no. | Cutting width |  |  | € |
|----------|---------------|---|---|---|
| GT412    | 40 mm         | 1/8   | 180   |   |

| Ref. no. | Description                                |  | € |
|----------|--|---|---|
| GT606    | Spare blade set for Akkurata/Vega          | 1/10/120  |   |
| GT606D   | Spare blade set for Akkurata/Vega with DLC | 1/10/120  |   |



All prices provided are the manufacturer's recommended retail prices and are inclusive of VAT. All information is subject to typing and printing errors. All prices are subject to change.

GT412=PI\_Vega\_II\_202412\_EN