

# Smart Collar™

A comfortable, well tolerated Elizabethan collar designed to protect cats and dogs from self-trauma during their recovery period.

- ✓ **Well tolerated**
- ✓ **Simple to fit**
- ✓ **Translucent**

An effective barrier to licking at sutures or self-trauma of wounds, the lightweight Smart Collar is made of flexible (but robust) translucent plastic with neoprene trim on the outer and inner edges. This protects furniture and owners' legs whilst allowing the free movement of the patients' head and neck so that they can walk, eat, drink and sleep in comfort aiding a quick recovery.

Available in 6 sizes, the Smart Collar is a practical solution to avoid self-trauma in cats and dogs of all breeds.



## YOUR OPTIONS

<b>A148-000</b>	Size 0 cat - Smart Collar Neck 14cm-22cm / Length 8cm
<b>A148-001</b>	Size 1 small - Smart Collar Neck 18cm-24cm / Length 9cm
<b>A148-002</b>	Size 2 medium - Smart Collar Neck 23cm - 30cm / Length 11cm
<b>A148-003</b>	Size 3 large - Smart Collar Neck 29cm-38cm / Length 16cm
<b>A148-004</b>	Size 4 x-large - Smart Collar Neck 35cm-52cm / Length 20cm
<b>A148-005</b>	Size 5 xx-large - Smart Collar Neck size 45cm - 64cm / Length 28cm



## How to get the right size - Smart Collar™ & SoftE Smart Collar™

Easily adjustable, measure animal's neck for correct size.

