


POLYPROPYLENE

BLUE



Type	Monofilament
Composition	Polypropylene, a polymer of propylene
Coating	None
Colour	Blue
Tissue reaction	Minimal
Absorption	Non absorbable
Presentation	Box with 12 sutures

*For quality
&
safety*

smi

www.sutures.be

POLYPROPYLENE

BLUE

CHARACTERISTICS

- ✓ High tensile strength, non absorbable
- ✓ Resistance to repeated bending
- ✓ Synthetic
- ✓ Monofilament with microscopic uniform diameter
- ✓ Pliability and softness
- ✓ Hermetically sealed packing

BENEFITS

Excellent and permanent tissue support

Use in cardiovascular surgery (heart prosthesis or vascular anastomosis)

Minimal tissue reaction

Soft passage through the tissues without sawing, tissue drag and trauma

Absence of capillarity (no support of bacterial growth)

Excellent knot tie down

Good handling properties

Easy and secure knotting

Guaranteed seal and product sterility

INDICATIONS General soft tissue closing
and/or ligation

especially in Cardiovascular surgery

Neurosurgery

Ophthalmic surgery

Microsurgery

Plastic surgery

Skin closure

Orthopaedics

Hernias (also as a mesh)

Gastrointestinal surgery

General surgery

Gynaecology

Obstetrics

Round bodied taper point




1/2 Circle



			Metric	USP	Length		Order number
	HR 10 mm		0.5	7/0	45 cm	blue	5050110
			0.7	6/0	75 cm	blue	5070110
	HR 15 mm		0.7	6/0	75 cm	blue	5070115
			1	5/0	75 cm	blue	5100115
			1.5	4/0	75 cm	blue	5150115
	HR 17 mm		1	5/0	75 cm	blue	5100117
			2	3/0	75 cm	blue	5200117
	HR 20 mm		1	5/0	75 cm	blue	5100120
			1.5	4/0	75 cm	blue	5150120
	HR 25 mm		2	3/0	75 cm	blue	5200120
			1	5/0	75 cm	blue	5100125
			1.5	4/0	75 cm	blue	5150125
			2	3/0	75 cm	blue	5200125
			3	2/0	75 cm	blue	5300125
			3.5	0	75 cm	blue	5350125
	HR 30 mm		4	1	75 cm	blue	5400125
			2	3/0	75 cm	blue	5200130
			3	2/0	75 cm	blue	5300130
			3.5	0	75 cm	blue	5350130
	HR 35 mm		4	1	75 cm	blue	5400130
			3	2/0	75 cm	blue	5300135
			3.5	0	75 cm	blue	5350135
	HR 40 mm		4	1	75 cm	blue	5400135
			3.5	0	75 cm	blue	5350140
			5	2	75 cm	blue	5500140

Round bodied taper point double armed




1/2 Circle




			Metric	USP	Length		Order number
	2 x HR 8 mm		0.4	8/0	60 cm	blue	5040108-2
			0.5	7/0	60 cm	blue	5050108-2
			0.7	6/0	60 cm	blue	5070108-2
	2 x HR 10 mm		0.5	7/0	60 cm	blue	5050110-2
			0.7	6/0	60 cm	blue	5070110-2
			1	5/0	60 cm	blue	5100110-2

Round bodied taper point double armed



			Metric	USP	Length		Order number
1/2 Circle 	2 x HR 13 mm		0.5	7/0	60 cm	blue	5050113-2
			0.7	6/0	60 cm	blue	5070113-2
			1	5/0	75 cm	blue	5100113-2
	2 x HR 15 mm		1.5	4/0	90 cm	blue	5150115-2
	2 x HR 16 mm		1	5/0	90 cm	blue	5100116-2
			1.5	4/0	90 cm	blue	5150116-2
			2	3/0	90 cm	blue	5200116-2
	2 x HR 20 mm		1	5/0	90 cm	blue	5100120-2
			1.5	4/0	90 cm	blue	5150120-2
			2	3/0	90 cm	blue	5200120-2
	2 x HR 25 mm		1.5	4/0	90 cm	blue	5150125-2
			2	3/0	90 cm	blue	5200125-2
			3	2/0	90 cm	blue	5300125-2
	2 x HR 30 mm		2	3/0	90 cm	blue	5200130-2
			3	2/0	90 cm	blue	5300130-2
	2 x HR 35 mm		3	2/0	90 cm	blue	5300135-2
	HR 15 mm HR 30 mm		2	3/0	90 cm	blue	5200115-30

Round bodied taper point double armed



			Metric	USP	Length		Order number
3/8 Circle 	2 x DR 8 mm		0.4	8/0	45 cm	blue	5040208-2
			0.5	7/0	60 cm	blue	5050208-2
	2 x DR 10 mm		0.7	6/0	60 cm	blue	5070210-2
	2 x DR 13 mm		0.7	6/0	60 cm	blue	5070213-2
	2 x DR 18 mm		1	5/0	90 cm	blue	5100218-2
			1.5	4/0	90 cm	blue	5150218-2

Round bodied taper point loop



Metric

USP

Length

Order number

1/2 Circle



HR 40 mm



3

2/0

150 cm

blue

5300940

Reverse cutting



Metric

USP

Length

Order number

3/8 Circle



DS 10 mm



0.7

6/0

45 cm

blue

5071510

DS 12 mm



0.5

7/0

45 cm

blue

5051512

0.7

6/0

45 cm

blue

5071512

1

5/0

45 cm

blue

5101512

1.5

4/0

75 cm

blue

5151512

DS 12 mm
Easy Pass



0.7

6/0

50 cm

blue

5071512P

1

5/0

50 cm

blue

5101512P

DS 16 mm



0.7

6/0

45 cm

blue

5071516

1

5/0

45 cm

blue

5101516

1.5

4/0

75 cm

blue

5151516

2

3/0

75 cm

blue

5201516

DS 16 mm
Easy Pass



0.7

6/0

50 cm

blue

5071516P

1

5/0

50 cm

blue

5101516P

1.5

4/0

50 cm

blue

5151516P

DS 19 mm



0.7

6/0

75 cm

blue

5071519

1

5/0

75 cm

blue

5101519

1.5

4/0

45 cm

blue

5151519-45

1.5

4/0

75 cm

blue

5151519

2

3/0

45 cm

blue

5201519-45

2

3/0

75 cm

blue

5201519

3

2/0

75 cm

blue

5301519

3.5

0

75 cm

blue

5351519

DS 19 mm
Easy Pass



1

5/0

50 cm

blue

5101519P

1.5

4/0

50 cm

blue

5151519P

2

3/0

50 cm

blue

5201519P

3

2/0

50 cm

blue

5301519P

Reverse cutting



3/8 Circle



		Metric	USP	Length		Order number
3/8 Circle	DS 24 mm	1.5	4/0	75 cm	blue	5151524
		2	3/0	45 cm	blue	5201524-45
		2	3/0	75 cm	blue	5201524
		3	2/0	45 cm	blue	5301524-45
		3	2/0	75 cm	blue	5301524
		3.5	0	75 cm	blue	5351524
		4	1	75 cm	blue	5401524
	DS 24 mm Easy Pass	2	3/0	50 cm	blue	5201524P
		3	2/0	50 cm	blue	5301524P
	DS 26 mm	3	2/0	45 cm	blue	5301526-45
		3	2/0	75 cm	blue	5301526
	DS 30 mm	2	3/0	75 cm	blue	5201530
		3	2/0	75 cm	blue	5301530
		3.5	0	75 cm	blue	5351530
		4	1	75 cm	blue	5401530
	DS 36 mm	2	3/0	75 cm	blue	5201536
		3	2/0	75 cm	blue	5301536
		3.5	0	75 cm	blue	5351536
		4	1	75 cm	blue	5401536
	DS 40 mm	2	3/0	75 cm	blue	5201540
		3	2/0	75 cm	blue	5301540
		3.5	0	75 cm	blue	5351540
		4	1	90 cm	blue	5401540
		5	2	90 cm	blue	5501540
	DS 45 mm	4	1	90 cm	blue	5401545
	DS 90 mm	4	1	90 cm	blue	5401590

Reverse cutting



1/2 Circle



		Metric	USP	Length		Order number
1/2 Circle	HS 22 mm	2	3/0	45 cm	blue	5201422
	HS 24 mm	3	2/0	75 cm	blue	5301424

Cutting 			Metric	USP	Length		Order number
Straight 	GS 40 mm 		3	2/0	75 cm	blue	5302140
			1.5	4/0	90 cm	blue	5152160
	GS 60 mm 		2	3/0	75 cm	blue	5202160
			3	2/0	75 cm	blue	5302160
	GS 75 mm 		3	2/0	75 cm	blue	5302175
Round bodied tapercut point double armed 			Metric	USP	Length		Order number
3/8 Circle 	2 x DRT 12 mm 		0.7	6/0	75 cm	blue	5070412-2
Round bodied tapercut point 			Metric	USP	Length		Order number
1/2 Circle 	HRT 10 mm 		0.7	6/0	75 cm	blue	5070510
	HRT 25 mm 		3	2/0	75 cm	blue	5300525
Round bodied tapercut point double armed 			Metric	USP	Length		Order number
1/2 Circle 	2 x HRT 16 mm 		1	5/0	75 cm	blue	5100516-2
	2 x HRT 22 mm 		1.5	4/0	90 cm	blue	5150522-2
	2 x HRT 25 mm 		1.5	4/0	90 cm	blue	5150525-2
			2	3/0	90 cm	blue	5200525-2
			3	2/0	90 cm	blue	5300525-2
	2 x HRT 30 mm 		2	3/0	90 cm	blue	5200530-2
			3	2/0	90 cm	blue	5300530-2

Soft cutting point coronary double armed



3/8 Circle



2 x CC 8 mm



Metric

USP

Length

Order number

0.4

8/0

60 cm

blue

5042608-2

0.5

7/0

75 cm

blue

5052608-2

0.7

6/0

75 cm

blue

5072608-2

2 x CC 10 m



0.5

7/0

75 cm

blue

5052610-2

0.7

6/0

75 cm

blue

5072610-2

2 x CC 13 mm



0.5

7/0

75 cm

blue

5052613-2

0.7

6/0

75 cm

blue

5072613-2

1

5/0

75 cm

blue

5102613-2

Soft cutting point coronary black double armed



3/8 Circle



2 x CC 8 mm black



0.4

8/0

60 cm

blue

5042608-2VB

0.5

7/0

75 cm

blue

5052608-2VB

0.7

6/0

75 cm

blue

5072608-2VB

2 x CC 13 mm black



0.7

6/0

75 cm

blue

5072613-2VB

1

5/0

75 cm

blue

5102613-2VB